


Phase fraction calculation of microstructure using ImageJ

BY


HARIHARAN V S

STEP-1




1. Open ImageJ software
2. Select File → Open

STEP-2




1. Select the required file
2. Click Open

STEP-3




The selected image appears

STEP-4




1. Select Image → Duplicate
2. Another copy of the image is created and file name prompt appears

STEP-5




1. Enter the name of duplicate file (which is of your wish)
2. Click OK
3. The duplicate image appears

STEP-6



Click on the rectangular box icon




STEP-7


The area other than the bottom unwanted area is selected

This area should not be counted for analysis

STEP-8




STEP-9





Select Image→Type→8-bit

STEP-10




Select Image→Adjust→Threshold

STEP-10




Adjust the scale to colour the phase of interest

STEP-11




Adjust the scale such that only the phase of interest is coloured red and then click APPLY

STEP-12



Select Analyze → Analyze Particles

STEP-13



1.Ensure that “summarize” option is checked

2.Click OK

Summary

Slice	Count	Total Area	Average Size	%Area	Mean
6-2.tif	6550	145.965	0.022	75.415	255

The phase fraction of shaded region is obtained

CAUTION

- ❖ The above procedure should be carried out for a large number of images taken in the same magnification and the average value should be reported.
- ❖ The obtained phase fraction will vary for different threshold values. Hence, it is necessary to choose the optimum threshold value which shades only the region/phase of interest.
- ❖ The images chosen for analyses should be clear and free of scratches.
- ❖ If the phases don't have a sharp contrast difference, it will be difficult to get accurate results.